Carboxylic Acid-Terminated Silicon Quantum Dots

Silicon Quantum Dots dry powder or dispersed in organic solvent

Description

Acid-terminated silicon quantum dots are readily dispersible in wide variety of polar organic solvents including toluene, methanol, ethanol, and buffered water. These particles can be used in photovoltaic devices, light emitting diodes, sensors, and optical film applications.

Advantages Over Traditional QDs

• Free of toxic metals (e.g., Cd, Pb, In) or phosphines

Applied Quantum Materials Inc.

- Bright PL, tunable from 600 to 950 nm (visible to near-IR)
- Low self-absorption due to large Stokes-shift >400 meV
- Stable PL at elevated temperatures > 100°C and high humidity

Product Specifications

	Size	PL _{max}	Catalog No.
Particle Sizes Available	2 nm	665 ±20 nm	14-0202-S2
	3 nm	710 ± 20 nm	14-0202-S3
	4 nm	780 ± 20 nm	14-0202-S4
	6 nm	845 ± 20 nm	14-0202-S6
	8 nm	980 ± 20 nm	14-0202-S8
Material Composition	Silicon		
Forms	Orange/red powders or solutions		
Photoluminescence	λ _{em} 600 to 1000 nm		
FWHM	<120 nm		
PL Lifetime	>50 µs		
Functionalization Method		QY	Shelf Life
Ι		10-40% ± 5%	12 months ¹
Ш		30-60% ± 5%	3 months ¹
¹ Note: The shelf life is approximate and requires proper storage conditions. Materials improperly stored can oxidize and lose optical properties.			





Uses & Handling Recommendations

- Shipped as powders or in solution. 1 mL, 5 mL, and 20 mL solution in glass vials (bulk can be supplied upon request).
- Typical concentrations ~ 3 mg/mL.
- Exposure of acid terminated SiQDs to water will cause oxidation. Minimize exposure to water whenever possible. Water soluble SiQDs can be made upon request.
- Sonication can be used to help disperse the acid terminated SiQDs in desired organic solvent (usually ethanol, methanol, and other polar organic solvents).

<u>Contact us</u> for purchasing/customization options. AQM can tailor the surface chemistry to provide SiQDs suitable for specific applications.

Characterization Data




Figure 1. (A) Photoluminescent emissions; (B) Particle size histogram; (C) Bright field TEM image.

Applied Quantum Materials Inc. B211, 2011 – 94 St. Edmonton, AB, Canada, T6N 1H1

Solutions in Silicon

contact@aqmaterials.com Tel: 1(587) 635-2060 www.aqmaterials.com